

# BIOLUX P-III

## Interim 12-month results for Superficial Femoral Artery (SFA)

### Conclusions

- BIOLUX P-III is a real all-comers registry evaluating the safety and the performance of Passeo-18 Lux Drug-Coated Balloon (DCB) in daily practice
- At 12 months the results of the first 131 SFA subjects demonstrate safety and performance comparable to the results reported in competitor global registries<sup>1</sup>
  - 91.8 % Freedom from Major Adverse Events (MAE)<sup>2</sup>
  - 94.0 % Freedom from clinically-driven Target Lesion Revascularization (cd-TLR)<sup>3</sup>
  - 82.1 % Primary Patency (PP)<sup>4</sup>
- High risk population and complex lesions subgroup analysis will reveal further insights on Passeo-18 Lux benefits in daily routine practice

### Study Design

Prospective, international, multi-center, all-comers registry investigating safety and efficacy data on the Passeo-18 Lux DCB in a real world population of subjects with atherosclerotic disease of the infrainguinal arteries in at least 700 patients

#### Principal investigator

- Prof. G. Tepe, Klinikum Rosenheim, Germany

### Endpoints

#### Primary endpoint

- Clinical: Freedom from MAE at 6 months
- Performance: Freedom from cd-TLR at 12 months

#### Secondary clinical endpoints (selected)

- Technical success
- Procedural success
- Device success
- Amputation-free survival at 6, 12 and 24 months
- PP rate at 12 and 24 months
- Freedom from cd-TLR at 6 and 24 months post index procedure
- Freedom from MAE at 12 and 24 months
- Clinical success defined as an improvement of Rutherford Classification (RC) at 6, 12 and 24 months follow-up of one class or more
- Changes in Ankle Brachial Index (ABI) measurements at 6, 12 and 24 months follow-up
- Patient-reported outcomes assessment: pain score, walking impairment questionnaire at 6, 12 and 24 months compared to the pre-procedure score

### Dedicated Subgroups

Diabetes; renal insufficiency; TASC II C&D lesions; heavily calcified lesions; RC 3+; 65 yrs+; Below-The-Knee (BTK); popliteal lesions, occlusions



## Baseline Characteristics

Age, yrs (mean ± SD), [Min; Max]	69.2 ± 10.3 [67.4 ; 70.9]
Male (n, %)	82 (62.6 %)

Medical history	n = 131	%
Hypertension	110	84.0 %
Hyperlipidemia	90	68.7 %
Smoking	103	78.6 %
Current smokers	43	41.7 %
History of PAOD	80	61.1 %
Previous PVI / surgeries	78	59.5 %
Diabetes	54	41.2 %
Coronary Artery Disease (CAD)	51	38.9 %
Cerebrovascular disease	34	26.0 %

## Lesion Characteristics


	n = 145	%
Lesion length, mm (mean ± SD)	80.4 ± 70.7	
Reference vessel diameter, mm (mean ± SD)	5.0 ± 0.8	
Diameter stenosis (%)	86.2 ± 13.4	
De novo lesion (n, %)	73	50.3 %
Occlusion (n, %)	32	22.1 %
In-stent restenosis (n, %)	21	14.5 %
Restenosis (n, %)	19	13.1 %

## Procedural Details

	n = 145	%
Lesions treated with DCB only	124	85.5 %
Lesions treated with DCB and a stent	21	14.5 %


## Freedom from MAE – SFA<sup>2,\*</sup>

(adjudicated by an independent CEC)




## Freedom from cd-TLR – SFA<sup>3,\*\*</sup>

(adjudicated by an independent CEC)



MAE <sup>2,*</sup>	n = 10
Death (procedure or device related death within 30 days post index procedure)	1
cd-TLR	8
Target limb major amputation	1

## PP – SFA<sup>4,\*\*</sup>



<sup>1</sup> ILLUMENATE Global Registry, INPACT Global, Lutonix Global SFA Real World Registry

<sup>2</sup> Major Adverse Event: Composite of device and procedure related mortality through 30 days post index procedure, major target limb amputation and clinically-driven Target Lesion Revascularization (cd-TLR). MAE are adjudicated by an independent Clinical Events Committee (CEC)

<sup>3</sup> Clinically-driven Target Lesion Revascularization (cd-TLR) is any re-intervention performed for ≥50 % diameter stenosis (visual estimate) at the target lesion after documentation of recurrent clinical symptoms of the patient

<sup>4</sup> Primary Patency (PP) is freedom from >50 % restenosis in the target lesion as indicated by a duplex ultrasound peak systolic velocity ratio (PSVR) >2.5 or by visual assessment of an angiogram with no clinically driven reintervention

\* Data presented by Prof. M. Brodmann at CIRSE 2016

\*\* Data presented by Dr. K. Keirse at LINC 2017