



# BIOFLEX PEACE

## 12-month results of full cohort

### Conclusions

- Full cohort 12-month\* results show Primary Patency (PP)<sup>1</sup> of 73.6%\* and Freedom from clinically driven Target Lesion Revascularization (Fcd-TLR) of 96.2%\* in this all-comers registry, which are in line with published study data for the Pulsar stent in controlled trials
- Clinical success of 80.8% (Improvement in Rutherford Class)
- Sub-group analysis has shown that patients treated with BMS and DCB have numerically higher PP than patients treated with BMS alone – a “Paclitaxel effect”
- Mean stent oversizing of 0.8 mm suggests that in the real-world, interventionalists are applying “minimal oversizing” to reduce the COF to least possible. Next step is to choose a stent design that reduces COF further

### Study Design

Prospective, multi-center, all-comers registry investigating safety and efficacy of the 4F Pulsar-18 stent in real world population in subjects with atherosclerotic disease of the femoropopliteal arteries

### Principal Investigator

- Dr. M. Lichtenberg, Arnsberg, Germany

### Endpoints

#### Primary endpoint


- 6-month Major Adverse Event (MAE) rate
- 12-month\* Primary Patency<sup>1</sup>

#### Secondary clinical endpoints (selected)

- Primary Patency at 6 and 24 months
- Fcd-TLR at 6, 12\* and 24 months
- Clinical success
  - Improvement of  $\geq 1$  Rutherford Class
  - Improvement in ABI

### Patient Demographics

		n = 164
Hypertension		146 (89.0%)
Dyslipidemia		131 (79.9%)
Smoker (current or previous)		119 (72.6%)
Diabetes		52 (31.7%)
Renal insufficiency		32 (19.5%)
<b>Critical Limb Ischemia (CLI)<sup>3</sup></b>		<b>29 (18.7%)</b>
Rutherford	0	1 (0.6%)
(n = 155)	1	3 (1.9%)
	2	59 (38.1%)
	3	63 (40.6%)
	4	15 (9.7%)
	5	12 (7.7%)
	6	2 (1.3%)
Ankle brachial index (n = 129)		0.66



### Lesion Characteristics

		n = 178
Lesion length (mm) (n = 143)		119.2 ± 103.69
Mean ref. vessel diameter		5.03 mm
Mean implanted strut diameter		5.83 mm
TASC A lesion		46 (25.8%)
TASC B lesion		60 (33.7%)
TASC C lesion		35 (19.7%)
TASC D lesion		37 (20.8%)
<b>TASC C/D lesions</b>		<b>40.5%</b>
Calcification	0 None	29 (16.3%)
	1 Mild	71 (39.9%)
	2 Moderate	43 (24.2%)
	3 Severe	35 (19.7%)
	Moderate & Severe	78 (43.8%)

## 12-month\* Results

Population	Endpoint	Vessel	n (patients)	A.L.L. (mm) mean ± SD	%
Full cohort	Freedom cd-TLR	Femoropopliteal	125/130	119.1 ± 103.2	96.2%
Imaging cohort <sup>4</sup>	PP	Femoropopliteal	67/91	138.5 ± 111.7	73.6%

## 12-month Results in Perspective

Study	Device	A.L.L.	PP	Fcd-TLR
BIOFLEX PEACE*	Pulsar-18	11.9 cm	73.6% <sup>5</sup> (67/91)	96.2% (125/130)
4EVER	Pulsar-18	10.8 cm	73.4%	85.2%
PEACE	Pulsar-18	11.2 cm	79.5%	81.0%
RESILIENT	Lifestent (C.R. Bard)	6.2 cm	81.3%	87.3%
ZILVER FLEX (PTX study)	Zilver Flex (Cook Medical)	6.3 cm	73.0%	77.0%
ZILVER PTX	Zilver PTX DES (Cook Medical)	6.6 cm	83.1%	90.5%
DURABILITY II	EverFlex (Medtronic)	8.9 cm	77.2%	n/a
SUPERA	Supera (Abbott)	9.0 cm	84.7%	n/a
DURABILITY	EverFlex (Medtronic)	9.6 cm	72.2%	79.1%
ABSOLUTE	Absolute (Abbott)	10.1 cm	63.0%	n/a

BIOFLEX PEACE 12-month data supports 4EVER and PEACE data in similar lesion length with comparable data to 6F BMS devices and Zilver PTX DES.


## Treatment Regime Influences Outcome?

Lesions	n = 178
Stent only	113 [63.5%]
Stent + DCB	48 [27.0%]
Stent + other device (thrombectomy, atherectomy, scoring balloon etc.)	17 [9.6%]

	Population	Endpoint	Vessel	n (lesions)	A.L.L. (mm) mean ± SD	%
Stent	Full cohort	Freedom cd-TLR	Femoropopliteal	89/90	121.2 ± 113.1	98.9%
	Imaging cohort <sup>4</sup>	PP	Femoropopliteal	38/54	137.3 ± 112.7	<b>70.4%</b>
Stent + DCB	Full cohort	Freedom cd-TLR	Femoropopliteal	34/35	99.8 ± 80.6	97.1%
	Imaging cohort <sup>4</sup>	PP	Femoropopliteal	23/28	94.6 ± 84.7	<b>82.1%</b>


## Clinical Success – Change in Rutherford

80.8% of subject show improvement of ≥1 Rutherford Class after 1 year<sup>6</sup>



Reference: Lichtenberg M., presented at Charing Cross Symposium April 2017

## Clinical Success – Change in ABI



\* 395 day results

1 Defined as freedom from >50% restenosis with no clinically-driven reintervention

2 Not for all subjects all measures are available

3 Rutherford classification 4-6 [29 out of 155 subjects]

4 Duplex data available. Duplex imaging was performed upon need and investigators discretion

5 Imaging cohort

6 Data available for 110 subjects

BIOTRONIK AG  
Ackerstrasse 6  
8180 Bülach · Switzerland  
Tel +41 (0) 44 8645111  
Fax +41 (0) 44 8645005  
info.vi@biotronik.com  
www.biotronik.com

© BIOTRONIK AG

All rights reserved. Specifications are subject to modification, revision and improvement.

 **BIOTRONIK**  
excellence for life